

Certificate of Conformance

**RAPICIDE® Test Strip**

**Part Number:** ML02-0120

**Lot Number:** 034105

**Expiration Date:** 2021-04-02

I certify that the RAPICIDE test strips indicate a positive reaction with a solution mixed to a concentration of 2.5% glutaraldehyde.

I certify that the RAPICIDE test strips indicate a negative reaction with a solution mixed to a concentration of 1.3% glutaraldehyde.


The statistical analysis and approved test procedure is on file and available for review.



Quality Analyst

14May2019

Analysis Date



Quality Assurance

14May2019

Release Date

**Medivators Inc.**

14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)



**MEDIVATORS**  
A Cantel Medical Company

Document  
Number: **50092-591**

Rev.  
Level **E**

**TITLE**

**FORM, RAPICIDE GLUT TS**

REV	DESCRIPTION	Effective Date	ORIGINATOR	ECO
D	ADD PASS/FAIL		SPW	32460
E	Update form to include test solutions, change info entry to lists/editable fields, correct spelling error, include kit and vial lot fields.		PBS	32975

**Part Number** ML02-0120

**Description** RapicideTest Strips

**Lot number** Vial:034105 Kit:034105

**Tested by** Nick Votrobeck

**Date** 14May2019

Note: Test form to be scanned and sent to incoming inspectors.

**Vial 1**

Concentration	Replicate 1	Replicate 2	Replicate 3
2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

**Vial 2**

2.5 %	Pass	Pass	Pass
1.3 %	Pass	Pass	Pass

**Visual Inspection**

**Testing Method: Unaided Eye**

	Pass / Fail	Comment
Cleanliness	Pass	Clean
Seal of caps	Pass	Sealed
Number of test strips -0/+3	Pass	60 per vial
Adherence of test strip pad	Pass	Adhered

**Testing Solutions**

Concentration	Lot# / Exp
2.5 %	415730 / 14Jun2019
1.3 %	CL19134-1 / 14Jun2019

*Nick Votrobeck* 14May2019

CONFIDENTIAL AND PROPRIETARY INFORMATION OF MEDIVATORS INC. COPYING AND DISSEMINATION PROHIBITED, UNLESS AUTHORIZED BY MEDIVATORS. RETURN REQUIRED UPON DEMAND.