

**Certificate of Analysis**

**RAPICIDE<sup>®</sup> PA**

**Part A**


Product Code: **ML02-115**

Lot Number: **036419**

Expiration Date: **2021/03**

Components	Specifications	Analysis
H <sub>2</sub> O <sub>2</sub>	20.0 – 24.0% (w/w)	<b>23.7 %</b>
PAA	4.0 – 6.0 % (w/w)	<b>5.3 %</b>
Acetic Acid	8.0 – 10.0% (w/w)	<b>9.0 %</b>
pH, (as is)	0.5 – 1.1 (as is)	<b>0.8</b>

This lot has been reviewed and meets specified requirements.

  
\_\_\_\_\_  
Quality Analyst  
\_\_\_\_\_  
Analysis Date  
\_\_\_\_\_  
Quality Assurance  
\_\_\_\_\_  
Release Date

**Medivators Inc.**  
14605 28th Avenue North  
Minneapolis, MN 55447 USA  
Telephone: (763) 553-3300  
Fax (763) 553-3387  
[www.medivators.com](http://www.medivators.com)

RAPICIDE<sup>®</sup> is a registered trademark of Medivators Inc.