

Certificate of Analysis

RAPICIDE® PA

Part B

Product Code: **ML02-0115**
Lot Number: **038318**
Expiration Date: **2020/05**

| Components | Specification | Analysis |
|---------------|---------------------|--------------------|
| Sodium | 743.0 - 909.0 mEq/L | 794.0 mEq/L |
| Benzotriazole | 0.9 – 1.1% w/w | 0.9 %w/w |
| pH | 11.4 to 12.4 as is | 11.6 |


This lot has been reviewed and meets specified requirements.



Quality Analyst

26-NOV-18

Analysis Date



Quality Assurance

30-NOV-18

Release Date

Medivators Inc.
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387

RAPICIDE® is a registered trademark of Medivators Inc.