



## Quality Control Report and Certificate of Analysis

Formula: BC-1-25

Lot #: 03XX07

### Concentrate Specification

|                      |                               |
|----------------------|-------------------------------|
| Sodium               | <u>1187.7</u> mEq/L $\pm$ 2 % |
| Chloride             | <u>402.6</u> mEq/L $\pm$ 2 %  |
| Bicarbonate          | <u>785.1</u> mEq/L $\pm$ 2 %  |
| pH                   | <u>7.0 - 8.0</u> pH units     |
| Conductivity @ 25 °C | <u>67</u> $\pm$ 3 mS/cm       |

When mixed using the correct equipment and dilution with AAMI certified water BC-1-25,  
Lot # 03XX07, meets the above specifications.

Approved for Release

Quality Control

Elaine Troy

Date:

10/11/07