



MINNTECH[®]
A Cantel Medical Company

Certificate of Analysis

77001-000

Lot Number:

626832

Volume:

21633

Liters

Chemical Manufacture Date:

03/02/09

<u>Components</u>	<u>Specification</u>	<u>Analysis</u>
Hydrogen Peroxide	21.0% - 24.0%	23.8%
Peracetic Acid	4.5% - 6.0%	5.3%
Acetic Acid	8.0% - 10.0%	9.4%
pH	0.5 - 1.1	0.7

This lot has been reviewed and meets specified requirements.

Elaine Trogen

Quality Analyst

03/16/09

Analysis Date

MINNTECH CORPORATION
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387