



MEDIVATORS[®]
REPROCESSING SYSTEMS

Certificate of Analysis

RAPICIDE[®]


Lot Number: 627029

Expiration Date: 03-2010

Quality Control Report

Analyte / Property	Test Method	Specification	Analysis
Appearance	Visual	Pale Yellow to Green Liquid	PASS
Glutaraldehyde	20600-752	2.5% to 2.6% (w/w)	2.6 %
Refractive Index	20501-108	1.3425 to 1.3445	1.3435
pH (as is)	20002-154	6.3 to 6.4 (as is)	6.3
Specific Gravity (25°C)	20600-834	1.010 to 1.016	1.016


This lot has been reviewed and meets specified requirements.



Quality Analyst

03-04-09

Analysis Date



Quality Assurance

03/04/09

Release Date

MINNTECH CORPORATION
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387