



MINNTECH®
A Cantel Medical Company

Certificate of Analysis

MINNCARE®

Part Number:	78399-664 (78325-300)
Lot Number:	640037
Package Configuration:	55 gal. drums
Package Net Volume:	40 gal./drum
Date of Packaging:	11/2009
Expiration Date:	11/2010

Components	Specification	Analysis	
Hydrogen Peroxide	22.0 ± 2.0%	23.3	%
Peracetic Acid	5.0 ± 1.0%	5.5	%
Acetic Acid	9.0 ± 1.0%	9.0	%
pH	0.8 ± 0.3	0.9	
Specific Gravity	1.090 to 1.140	1.114	

Product has been filtered to $\leq 0.2 \mu\text{M}$ prior to final packaging.

This lot has been reviewed and meets specified requirements.

Quality Analyst

10-29-2009

Analysis Date

Quality Assurance

10-29-2009

Release Date

MINNTECH CORPORATION
14605 28th Avenue North
Plymouth, Minnesota 55447
Telephone: (763) 553-3300
Fax (763) 553-3387