



**MEDIVATORS**  
A Cantel Medical Company

## Certificate of Conformance

# RAPICIDE<sup>®</sup> PA

## Test Strips

**Part Number:**

**ML02-0129**

**Vial Lot Number:**

OG155065

**Expiration Date:**

07/2017

### Quality Control Report:

MEDIVATORS certifies that the test strips indicate a positive reaction with a solution mixed to a concentration of  $1000 \pm 30$  ppm Peracetic acid.

MEDIVATORS certifies that the test strips did not indicate a positive reaction with a solution mixed to a concentration of  $500 \pm 30$  ppm Peracetic acid.

This lot has been reviewed and meets specified requirements.



Quality Assurance

12Nov2015  
Release Date

Medivators Inc.  
14605 28<sup>th</sup> Avenue North  
Minneapolis, MN 55447 USA  
Phone: 1.800.444.4728  
Fax: 1.763.553.3387  
www.medivators.com