

*Certificate of Conformance*

**ACTRIL<sup>®</sup> Indicator Test Strips**

**Part Number:** **78280-000**

**Lot Number:** PA19086

**Expiration Date:** 2022-03-27

**Quality Control Report:**

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of  $50 \pm 10$  ppm Peracetic acid.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of  $5 \pm 1$  ppm Peracetic acid.

This lot has been reviewed and meets specified requirements.

  
Quality Assurance

4/17/2019  
Release Date