

Certificate of Conformance

ACTRIL[®] Indicator Test Strips

Part Number: **78280-000**

Lot Number: PA19091

Expiration Date: 2022-04-01

Quality Control Report:

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 50 ± 10 ppm Peracetic acid.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 5 ± 1 ppm Peracetic acid.

This lot has been reviewed and meets specified requirements.



Quality Assurance

4/29/2019

Release Date