

Certificate of Conformance

ACTRIL[®] Indicator Test Strips

Part Number: **78280-000**

Lot Number: PA19130


Expiration Date: 2022-05-10

Quality Control Report:

Minntech certifies that the test strips indicate a positive or pass reaction with a solution mixed to a concentration of 50 ± 10 ppm Peracetic acid.

Minntech certifies that the test strips indicate a negative or fail reaction with a solution mixed to a concentration of 5 ± 1 ppm Peracetic acid.

This lot has been reviewed and meets specified requirements.


Quality Assurance

6/13/2019
Release Date