

SCOPE BUDDY® COMPETENCY FORM



Facility Name _____

City / State _____

Contact Name _____

Telephone: _____

Employee will demonstrate competency in using the Medivators Scope Buddy

#	Demonstrated Skills	Pass	Redo	Comments
1	Scope Buddy setup and familiarization			
2	Tubing segment identification			
3	System operation (Start/Stop, Time Set arrow keys)			
4	Use of Scope Buddy Rinse Matrix for selection of correct endoscope flush time			
5	Appropriate choice of Scope Buddy hook-ups for the endoscope being flushed (DSD vs. universal tubings)			
6	Appropriate dilution for detergent used			
7	Endoscope flushing procedure using Scope Buddy			
8	Knowledge of special instructions for an ERCP endoscope elevator wire channel			
9	Daily Flow Validation Test			
10	Daily Decontamination Procedure			
11	Recommended overnight and weekend storage instructions			
12	Where / who to call for Medivators Customer Support			

Reviewer's Evaluation Results:

Pass _____

Redo _____

Employee's Name _____

Signature: _____

Date: _____

Title: _____

Reviewer's Name: _____

Signature: _____

Date: _____

Title: _____